

Child Growth Monitoring System : CGMS Series



SALIENT FEATURES :

- ◆ It gives the result of a child comprising of all three growth charts published by World Health Organisation (WHO), namely, UNDERWEIGHT, STUNTING and WASTING. It gives the CATEGORY of Underweight, Stunting and Wasting as well as shows the amount of Underweight and Stunting of a child.
- ◆ It can store the data of 500 children.
- ◆ HOST USB Port is provided by which the data of children can be copied on a PEN DRIVE and the same can be transferred to computer so that authentic malnutrition status can be provided at desired level of administration in a short time frame.

MAIN FEATURES :

- ◆ Computer software, "MAA" (Malnutrition Assessment and Analysis) by which one can analyse the various data about Growth Monitoring in tabular and Graphical form. One can view information of an individual child as well as Anganwadi, Place, Project, Block, District and State, individually or collectively.
- ◆ It can be equipped with rechargeable battery back-up by which the system can be operated for minimum 8 hrs on battery.
- ◆ It is equipped with very sturdy and large size of ABS bowl (As per Weights & Measures Rules) which can accommodate children ageing from 0 to 72 months.
- ◆ 25-key feather touch membrane keyboard with audible indication.
- ◆ 18-digit easy-to-read bright red LED display.
- ◆ Equipped with Real Time Clock, provides Date of each weighment.
- ◆ User friendly local language keyboard, keeping in view utilisation by Anganwadi workers. ◆ Portable and sturdy design.

ACCESSORIES :

- ◆ Rechargeable 6V, 4AH, SMF battery with battery charging circuit.
- ◆ "MAA" (Malnutrition Assessment and Analysis) Software.

Screen Shot of 'MAA' Software

